SR 200 APR Full-face Mask

Chemical Cartridge Application
PAPR Application

User's Instruction • Manuel de l'utilisation • Instrucciones de uso

Sundstrom Safety Inc.

143 Maple St. Warwick, RI 02888

Office: 1-401-434-7300 Toll Free: 1-877-SUNDSTROM Fax: 1-401-434-8300

info@srsafety.com www.srsafety.com

English	Please read and save these instructions	3
	Illustrations	29
Français	Prière de lire et de conserver	11
	Figures	29
Español	Lea y conserve estas instrucciones por favor ...	20
	Ilustraciones	29

Table of contents

- I NIOSH Approval Limitations
- II General Respirator Use Limitations
- III Introduction
- IV Configuration and selection guide
- V User Training and Instructions
 - 1. Changing the Protective Elements
 - 2. Inspection Prior to Use
 - 3. Donning/Fit Check/DoFFing
- VI Cleaning and Storage
- VII Spare parts
- VIII Warranty
- IX Illustrations
- X NIOSH/DHHS Labels (Separate insert pages included with this booklet)

I. NIOSH Approval Limitations

Note: All cautions and limitations do not apply to all applications. Refer to the NIOSH matrix approval label insert (add part number) to verify the applicable cautions and limitations.

Chemical cartridge APR application:

- A- Not for use in atmospheres containing less than 19.5 percent oxygen.
- B- Not for use in atmospheres immediately dangerous to life or health.
- C- Do not exceed maximum use concentrations established by regulatory standards.
- H- Follow established cartridge and canister change out schedules or observe ESLI to ensure that cartridges and canisters are replaced before breakthrough occurs.
- I- Contains electrical parts that may cause an ignition source in flammable or explosive atmospheres.
- J- Failure to properly use and maintain this product could result in injury or death.
- L- Follow the manufacturer's User's Instructions for changing canisters.
- M- All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.
- N- Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O- Refer to User's Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P- NIOSH does not evaluate respirators for use as surgical masks.

FF- Respirators are to be fit tested prior to use with the heaviest cartridges, canisters, filters and/or accessories intended to be used. Fit testing should also be conducted while wearing all personal protective equipment intended to be used. See User's Instructions for fit test requirements.

S- Special or critical User's Instructions and/or specific use limitations apply.

The following are NIOSH Cautions and Limitations for PAPP applications:

- A- Not for use in atmospheres containing less than 19.5 percent oxygen.
- B- Not for use in atmospheres immediately dangerous to life or health.
- C- Do not exceed maximum use concentrations established by regulatory standards.
- F- Do not use powered air-purifying respirators if air flow is less than four cfm (115 lpm) for tight fitting facepieces or six cfm (170 lpm) for hoods and/or helmets.
- H- Follow established cartridge and cartridge change schedules or observe ESLI to ensure cartridges and canisters are replaced before breakthrough occurs.
- I- Contains electrical parts that may cause an ignition source in flammable or explosive atmospheres.
- J- Failure to properly use and maintain this product could result in injury or death.
- L- Follow the manufacturer's User's Instructions for changing canisters.
- M- All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.
- N- Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O- Refer to User's Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P- NIOSH does not evaluate respirators for use as surgical masks.
- FF- Respirators are to be fit tested prior to use with the heaviest cartridges, canisters, filters and/or accessories intended to be used. Fit testing should also be conducted while wearing all personal protective equipment intended to be used. See User's Instructions for fit test requirements.
- S- Special or critical User's Instructions and/or specific use limitations apply.

II. General respirator use limitations

Full-face Air-Purifying Respirators (APR) are designed, tested and government approved to provide respiratory protection to a properly trained person in certain hazardous atmospheres.

The wearer must not enter a hazardous atmosphere that meets or exceeds the following limitations:

Maximum Use Concentrations (MUC)

- a. 50 times the exposure limit (TWA) for the contaminant present.
- b. Immediately dangerous to life or health (IDLH) concentrations for the contaminant.
- c. The limitations outlined in the applicable NIOSH approval.
- d. Any applicable limitations contained in the standards established by a regulatory agency (such as federal, state and local agencies) with jurisdiction over the wearer.

WARNING

This respirator cannot remove all the contaminants from a hazardous atmosphere. Pay close attention to all Maximum Use Concentrations (MUC), Limitations and other restrictions.

In addition to these Maximum Use Concentrations (MUC), Limitations, and the warnings and cautions placed in the text, the following general safety points must be read and understood by everyone intending to use the Sundström respirator.

- Do not use for paints containing Isocyanates.
- Do not use for protection against fumigants.
- This respirator is designed to be used only by TRAINED QUALIFIED PERSONNEL.
- Respirator selection shall be based upon full knowledge of, and professional evaluations of, the contaminants and their maximum content in the subject atmosphere.
- Air-purifying respirators DO NOT supply oxygen and MUST NOT BE USED in atmospheres containing less than 19.5 % oxygen by volume.
- This respirator MUST NOT BE USED for protection against any substance that may cause an allergic reaction, that contain Isocyanates or other restricted contaminants or that lack adequate warning properties.
- Air-purifying respirators shall not be used for structural fire fighting.

- DO NOT USE this respirator for protection in any atmospheres that are, or may become immediately Dangerous to Life or Health (IDLH), or in conditions that pose an immediate threat by air contaminants that are likely to have adverse delayed effects on life or health.
- This respirator MUST NOT BE USED for protection against substances that lack adequate warning properties.
- This respirator is not approved for protection from abrasive blasting, including sandblasting operations.
- This respirator must be equipped with the appropriate NIOSH approved air-purifying filters, chemical cartridges, or filter/cartridge combinations that are properly installed on the respirator.
- This respirator must be properly fitted to the individual who will wear it. Do not use the Sundström respirator if you have any hair growth between the skin and facepiece sealing surface such as stubble, beard growth, beard, moustache, or sideburns which cross the respirator surface (see OSHA standard 29 CFR 1910.134 Appendix A.I .A.9 for further clarification). Do not use this respirator if scars or other physical characteristics may interfere with a proper fit of the respirator.
- Always read cartridge and filter approval labels prior to use. Air-purifying elements designed only for protection against airborne particles MUST NOT BE USED for protection against vapors or gases. Air-purifying elements designed only for protection against vapors or gases MUST NOT BE USED for protection against airborne particulates.
- The useful service life of the air-purifying elements will vary with the concentration and nature of the contaminant(s) and the activity of the respirator wearer. Replace the air-purifying elements when you detect any indication of damage or loss of effectiveness, or when instructed by competent authority.
- If you experience any one of the following danger signals leave the hazardous area immediately! Return to fresh air, and seek assistance.
 - a. Breathing difficulty.
 - b. Dizziness or nausea.
 - c. You smell or taste the contaminants.
 - d. Any other noticeable physical problem.

If the appropriate protection is not specified by a competent authority, or should you feel uncertain about the proper application, selection, or maintenance of this equipment, consult your supervisor, safety professional, or local occupational health organization. Do not enter a suspected hazardous area. For answers to your questions contact Sundstrom Safety Inc. at 1-877-SUNDSTROM.

III. Introduction

The Sundström SR 200 Full-face Air-purifying Respirator (APR) has been tested and approved by NIOSH to provide the approved level of respiratory protection to a user.

Provides limited eye protection (ANSI Z87.1-2010).

The respirator consists of an outer mask with polycarbonate visor that covers the user's face, an integrated inner mask with inhalation and exhalation valves that covers the user's nose, mouth and chin, a rubber strap head harness with 6 mounting points that hold the respirator in place, and a filter adapter for connecting standard Sundström protective elements.

The inhaled air flows through a protective element and inhalation membrane into the inner mask. Part of the air flows past the inside of the visor in order to prevent misting. The exhaled air is discharged from the face piece through two exhalation valves. A laminated glass visor, a test adapter and spectacle frames for fitting corrective lenses are available as accessories.

This manual is divided into ten sections of equal importance. It will provide a respirator user with the basic information to properly use, inspect, don, doff, and perform maintenance for the Sundström respirator.

This manual provides only the basic elements of a complete respirator program. A well-designed respirator program also includes fit testing, proper documentation and active management by your local safety professional. All regulations and recommendations found in OSHA documents, especially 29 CFR 1910.134; EPA Standards; American National Standards Institutes (ANSI) Standard Z88.2; and all applicable state, provincial and local regulations shall be referenced in the complete safety program.

Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

This respirator should not be used without thorough understanding of all the components and their purpose. Read this entire booklet carefully and refer to it whenever a question arises. You must also satisfactorily complete all regulatory and local training requirements before entering a hazardous atmosphere.

Failure to read, understand and carefully follow these instructions prior to using this product may result in personal injury or death.

IV. Configuration and selection guide

FULL-FACE RESPIRATOR

Configuration, fig. 1

Component	Order Number
1. Polycarbonate visor	R01-1201
2. Upper frame half with screws	R01-1202
3. Head harness assembly	R01-1203
4. Set of membranes	R01-1204
a) Exhalation membranes, two	
b) Valve covers, two	
c) Inhalation membranes, three	
d) Dowels, two	
5. Pre-filter holder	R01-0605
6. Filter adapter SR 280-3	H09-0212
7. Gasket for filter connection	R01-1205
8. Carry Strap	R01-1206

FILTER/CARTRIDGES/PRE-FILTER Order Number

OV Cartridge, SR 218-6	H02-2121
CL/HC/SD/FM Cartridge, SR 231	H02-3621
OV/SD/CL/HC/HF Cartridge, SR 232	H02-3721
AM/MA Cartridge, SR 229	H02-4321
P100 filter SR 510	H02-1321
Pre-filter, SR 221	H02-0121

ACCESSORIES

	Order Number
Cleaning Wipes SR 5226	H09-0401
Test adapter	T01-1202
Permanent ID tag	R09-0101
Visor, laminated glass	T01-1203
Spectacle kit for corrective lenses SR 341	T01-1201
Decontamination Cap	R03-1406
Spark Arrester, SR 336	T01-2001
Welding Cassette, SR 84	T01-1212
Protective Film, SR 343	T01-1204
Protective Film, SR 353	T01-1205
Rubber Head Harness	T01-1215
Protective Hood, SR 64	H09-0321
Protective Hood, SR 345	H09-1021
SmallTalk Voice Amplifier, ST2-SR	T01-1219

Selection of protective elements

The protective elements approved for use on the Sundström Full-face Respirator are (4) chemical cartridges and (1) particulate filter. A pre-filter is also approved for optional use with all protective elements.

You can identify various protective elements by the color and protection designation of the NIOSH filter/cartridge label.

MAGENTA	P100
BLACK	OV
YELLOW	OV / SD / CL / HC / HF
WHITE	CL / HC / SD / FM
GREEN	AM / MA

P100 – Particulate filter (99.97 % filter efficiency level) effective against all particulate aerosols.

- OV** – Organic vapor
- SD** – Sulfur dioxide
- CL** – Chlorine
- HC** – Hydrogen chloride
- HF** – Hydrogen fluoride
- FM** – Formaldehyde
- AM** – Ammonia
- MA** – Methylamine

Chemical Cartridges

Each chemical cartridge is designed to provide respiratory protection against specific contaminants. A chemical cartridge absorbs and/or adsorbs specific vapors and gases from a contaminated atmosphere. This process continues until the sorbent becomes saturated and allows the contaminant to break through. OSHA requires a chemical cartridge change out schedule based on workplace aerosol concentration and no longer relies on breakthrough sensory detection to initiate chemical cartridge changing.

Particulate filters

Filters provide protection against particulates. All Sundström particulate filters trap and hold particles in the filtering media. As the amount of the captured contaminant in the media increases, breathing resistance also increases. Replace the filter when the change in breathing resistance becomes noticeable.

The P100 filter is 99.97 % efficient against all particulate aerosols.

Sundström Safety has conducted degradation testing of our P100 filter. Our test results show that even after an extensive period exposed to oil based contaminants this filter will not degrade.

Combination Protective Elements

An appropriate chemical cartridge and P100 filter can be combined to provide protection from a combination of certain vapors, gases and particulates.

Pre-filter

The Sundström pre-filter is not an protective element and can never be used as primary protection or as a substitute for the P100.

It is designed to prevent nuisance particulates from reaching the protective elements. This increases the life span of the primary protective element.

FOR USE WITH SR 500 POWERED AIR-PURIFYING RESPIRATOR (PAPR)

SR 200 Full-face respirator are approved for use with SR 500 PAPR fan unit assembly. Refer to the NIOSH approval label for approved configurations P/NL67-0421. In addition to this manual, the SR 500 manual P/N L67-0121 must be carefully read and followed by all persons who have, or will have, the responsibility for using or servicing the respirator.

FOR USE WITH SR 700 POWERED AIR-PURIFYING RESPIRATOR (PAPR)

SR 200 Full-face respirator are approved for use with SR 700 PAPR fan unit assembly. Refer to the NIOSH approval label for approved configurations P/NL67-0621. In addition to this manual, the SR 700 manual P/N L67-0721 must be carefully read and followed by all persons who have, or will have, the responsibility for using or servicing the respirator.

User Training and Instructions

WARNING

Do not complete the user training without reading and understanding all instructions in this booklet.

This training section is designed for the user and their safety manager. It contains the minimum information required to safely instruct users on how to don, doff, test, and change protective elements for the Sundström respirator. This training procedure is included as a supplement to the instructions contained in the remainder of this manual.

Local conditions may require additional instructions. Sundström encourages the safety professional responsible for the respiratory program to add all local information possible to increase worker safety and to supplement this basic text. You can obtain an Instructors Guide from your nearest Sundström Distributor.

1. WHY USE RESPIRATORY PROTECTION?

A. Why is breathing protection necessary?

Breathing protection is necessary whenever you are in an atmosphere that could injure your health. Depending on the substance present in the contaminated atmosphere you could be at risk of:

LUNG DAMAGE
Bronchitis, edema, pneumonia

DAMAGE TO VITAL ORGANS
Kidneys, liver, brain

CANCER

PERMANENT DISABILITY
asthma, silicosis, asbestosis or emphysema

or at worst
DEATH

WARNING

Do not enter a hazardous or unknown atmosphere without asking your safety professional what type of protection is required. Use all the protective equipment and wear it properly. At the first sign of trouble exit the hazardous area, seek help and report the problem to the proper authority.

B. What is a full face respirator?

The Sundström full face respirator consists of an outer mask with polycarbonate or glass visor that covers your face, an integrated inner mask with inhalation and exhalation valves that covers your nose, mouth and chin, a rubber strap head harness with 6 mounting points that hold the respirator in place, and a filter adapter for connecting standard Sundström protective elements.

C. What type of protection will a full face respirator provide?

A properly fitted full-face respirator can protect you against specific levels of some:

Particulates
Gases and Vapors
Viruses and Bacteria
Mold and Fiber

BUT full-face respirators **SHOULD NOT** be used:
In oxygen deficient areas

In atmospheres Immediately Dangerous to Life and Health (IDLH)

When the contaminants are unknown.

WARNING

Know your Maximum Use Limitations and follow the instructions of your safety professional.

D. How do I know which filter or cartridge to use?

Your safety professional is responsible for providing you information concerning what filter, cartridge or combination will provide adequate protection for the hazards at your worksite. You can identify various 10 protective elements by the color of the NIOSH label.

MAGENTA

Particles, Dust, Mists and Fumes

BLACK

Organic Vapors

YELLOW

Organic Vapors, Sulfur Dioxide, Chlorine, Hydrogen Chloride, Hydrogen Fluoride

WHITE

Chlorine, Hydrogen Chloride, Sulfur Dioxide, Formaldehyde

GREEN

Ammonia and Methylamine

The NIOSH/DHHS approval label for each product lists limitations and requirements. All current NIOSH/DHHS labels can be found in the respirator box printed on a sheet titled Section "X NIOSH/DHHS LABELS" of this manual. The Sundström P100 filter can be combined with each chemical cartridge to give the user combined vapor, gas and particulate protection.

E. What does the pre-filter do?

The Sundström pre-filter can only be used with a cartridge or P100 filter to reduce nuisance dusts. This extends the life of the protective element by capturing these large particulates, preventing premature loading of the primary protective elements. This pre-filter is not a protective element.

WARNING

The Sundström pre-filter is a thin filter that does not have any printing on either side. Do not use this filter as a protective element.

F. What if I can't smell anything, or I don't mind the smell?

The physiological damage from some contaminants can be long term and many take months to appear. The danger from these hazards may be masked by an unfamiliar or even pleasant odor. A certain smell or odor may not bother you, or cause you discomfort, but the contaminant can still be hazardous.

G. When do I change the filter or cartridge?

Chemical Cartridge – Loss of protection is frequently detected when the user experiences BREAKTHROUGH.

If you experience any one of the following danger signals leave the hazardous area immediately, return to fresh air. Check your respirator, and replace the protective element if necessary; and seek any necessary assistance.

- Breathing difficulty, dizziness or nausea, smell
- Dizziness or nausea

- c. You smell or taste the contaminants
- d. Any other noticeable physical effect.

If you smell, feel or taste anything unusual leave the contaminated area, check your respirator, and replace the protective element if necessary. In some atmospheres breakthrough may be hard to detect or may occur too late to provide you a safe exit from the contaminated atmosphere. An air-purifying respirator is not appropriate for use. However under certain conditions, the chemical cartridges can be used but must be changed after a specific exposure time. This schedule is set by a competent authority and controlled by your safety professional.

WARNING

Warning. An exhausted chemical cartridge will allow hazardous fumes or vapors to pass through to the user who may not notice any change in breathing resistance. Follow local procedures and all warnings contained in the users manual.

Filters – Filters clean the air of particles by capturing the solid particles and allowing breathable air to enter the respirator. As the filter loads up with these particulates inhalation becomes uncomfortable. When you detect increased inhalation resistance, leave the contaminated area and change the pre-filter, or if necessary, the protective filter.

Time use restrictions may also apply to some filters when efficiency declines or when they are exposed to specific contaminants. Filters must be replaced every time these limits are reached.

Pre-filter – This inexpensive nuisance pre-filter is used with a cartridge or P100 filter. When breathing resistance increases change the pre-filter first. If the resistance does not drop change the protective element.

1. Changing the protective elements

A. Installation / removal

Select the appropriate protection. Remain in a safe atmosphere and select the proper cartridge and appropriate filter.

B. Assembling the cartridge, P100 filter combination. Fig. 2.

- Place the P100 filter on top of the cartridge. Grasp both protective elements.
- Squeeze hard until you hear the P100 snap onto the cartridge.
- Pull on the filter and cartridge to check that they are properly connected.

- Place a pre-filter into the pre-filter holder.
- Snap the pre-filter holder onto the filter.
- Check to see that it is firmly attached to the cartridge or filter.

Installing the protective element into the respirator.

- Hold the protective element in your hand with the arrows on the label pointing towards the face. Insert one part of the large lip of the element into the filter adapter flange.
- Grab the flange with your finger and pull the lip of the adapter flange around the protective element until the entire protective element is seated inside the filter adapter.
- Squeeze the top edge of the element and the bottom edge of the respirator filter adapter together until the filter adapter flange is seated against the edge of the cartridge.
- Make sure that the gasket in the respirator connector is in place, screw the filter adapter with the protective element onto the respirator connector. Carefully tighten the filter adapter by hand to prevent damage to the thread and to ensure a good seal against the gasket.

C. Removing the protective elements from the respirator.

- Screw off the filter adapter with the protective element from the respirator connector.
- On the top of the respirator filter adapter flange is a small tab that touches the side of the protective element.
- Grasp the protective element with one hand and this tab with the finger and thumb of the other hand.
- Peel back on the tab while pulling on the cartridge.
- If you are using a stacked cartridge and filter, you must now remove the filter from the cartridge.
- Take a coin, or similar item, and place it in the space between the lower lip on the P100 Filter and the small raised plastic tab on the side of the cartridge. Twist the coin until the filter pops off the cartridge. Fig. 3.

NOTE

All cartridges and filters should be disposed of as hazardous waste in accordance with national waste disposal regulations.

2. Inspection prior to use

Before donning the respirator always perform the following checks.

- Ensure that the respirator is clean and dry.
- Check the respirator, valves, harness and visor for cracks, missing parts, and other defects.
- Check that the appropriate chemical cartridge, filter or combination is installed properly.

3. Donning/Fit check/Doffing

A. Donning

- Remove any hood, glasses or ear protection from head.
- Slacken the four elastic straps by moving the strap holders forward at the same time pulling the straps. Fig. 4. Slacken the upper two inelastic straps by opening the buckles.
- Move the head harness upwards, place your chin in the facepiece chin support and pull the head harness over your head. Fig. 5.
- Tension the elastic straps in pairs by pulling the free strap ends towards the rear. Fig. 6.
- Adjust the lengths of the upper pair of straps and fix by means of the buckles.
- Adjust the fit of the respirator on your face, so that it fits firmly but comfortably.
- Perform a fit check.

B. Fit Check

WARNING

A proper fit check must be successfully completed each time the respirator is worn. Failure to achieve a proper face-to-respirator seal before you enter a hazardous atmosphere may result in little or no respiratory protection.

- Install a protective element.
- Snap the pre-filter holder onto the protective element.
- Donning the respirator.
- Place the palm of your hand lightly over the hole on the pre-filter holder to make it tight. Fig. 20.
NOTE! Do not push so hard that the respirator's shape is affected.
- Inhale gently, the respirator should pull inward on the face.
- Hold your breath, the respirator should remain collapsed on the face.
- Remove the palm of your hand from the pre-filter holder.
- Breathe normally.

If any leakage is detected, check the inhalation and exhalation valves or adjust the straps of head harness. Repeat the fit check until there is no leakage.

C. Doffing

- Slacken the four elastic straps in pairs by moving the strap holder forward. The two inelastic straps need not be released. Fig. 7.
- Pull the head harness forward over your head and remove the respirator.
- Decontaminate, clean and store respirator as required.

SR 500 PAPR application

For donning and doffing Instructions, refer to the SR 500 PAPR User's Instructions P/N L67-0121.

SR 700 PAPR application

For donning and doffing Instructions, refer to the SR 700 PAPR User's Instructions P/N L67-0721.

VI Cleaning and storage

A. Disassembly

- Remove the filter adapter with protective element.
- Remove the covers from the exhalation valves and remove the membranes (two).
- Remove the inhalation membranes (three).
- Remove the head harness.
- If necessary remove the visor. See section VII, A.

B. Cleaning and assembling

A dirty respirator causes reduced protection. A written program for cleaning and sanitizing should be established and followed. This program should include all federal and local requirements.

- Clean with warm soapy water or other respirator cleaner and a soft brush, lint free cloth or sponge. Do not exceed 104 °F (40 °C). Never use solvents for cleaning.
- Sanitize the respirator. Use any approved sanitizing or disinfecting agent or a solution of 2 tablespoons of household bleach mixed in one gallon of water (50-PPM chlorine); or 1 teaspoon of tincture of iodine (50-PPM iodine) in one gallon of water. Immerse the respirator in the solution for at least 2 minutes.
- Rinse thoroughly.
- Check the condition of the respirator and all the parts. Replace any damaged or defective item.
- Shake the respirator to remove excess water and allow to dry.
- Perform a quick fit check.

C. Storage

Store the respirator, filters and cartridges in a clean, cool, dry place out of direct sunlight.

Shelf life for protective elements in unopened packages: Particulate filter 10 years, chemical cartridge 5 years.

The expiration date is printed on the label: month, year, e.g. May 2032.

VII Spare parts

Use only genuine Sundström parts. Don't modify the equipment. The use of "pirate parts" or any modifications may reduce the protective function and will compromise the approvals granted to the product.

A. To change the visor

The visor is mounted in a groove running around the visor opening of the outer mask and is held in place by one upper and one lower frame half.

- Use a 2.5 mm Allen key to remove the two screws holding the frame halves together. Fig. 8.
- Carefully remove the upper frame half. Fig. 9.
- Carefully pry the top part of the mask off the visor, and remove the visor from the lower groove. Fig. 10,11. Take this opportunity to clean the groove, if necessary.
- Markings are made to show the centres of the visor, frame halves and mask. Fig. 12. Press the new visor into the groove, making sure that the centre markings are in line. To make assembly easier, coat the slot with a soap solution or similar liquid.
- Carefully prise the top half of the facepiece over the visor, and make sure that the visor is in the groove in the facepiece. Fig. 10.
- Prise the upper frame half, making sure that the centre markings are in line. Fig. 12.
- Fit the screws and tighten them alternately until the two halves of the frame are firmly in contact. Fig. 8.
- Perform a fit check.

B. To change the inhalation membranes (three)

- One membrane is in the centre of the inner mask on a fixed dowel. Prise off the membrane and fit a new membrane. Fig. 13.
- Two membranes are fitted, i.e. one on each inside of the inner mask. The dowels for these membranes are removable and should be changed whenever the membranes are changed.
- Prise off the membranes and dowels.
- Prise the new membranes onto the dowels. Fig. 14.
- The membranes should rest on the larger flange, i.e. thread the dowel with the membrane from the inside of the facepiece, through the valve seat, with the smaller flange first. Fig. 15.

C. To change the exhalation membranes (two)

The exhalation membranes are mounted on a fixed dowel on the inside of the valve covers on each side of the outer facepiece. The covers should be changed whenever the membranes are changed.

- Snap the valve covers off the valve seats. Fig. 16.
- Prise off the membranes. Fig. 17.
- Press the new membranes onto the dowels. Carefully check that the membranes are in contact with the valve seats all round.
- Press the valve covers into place. A clicking sound indicates that the cover has snapped into place.
- Perform a fit check.

D. To change the head harness

The head harness can be ordered as a spare part only as a complete harness.

- Snap the strap holders of the head harness off the facepiece strap mountings. Fig. 18, 19.
- Check that the straps are not twisted and fit the new head harness.

VIII Warranty

Supplier warrants to each purchaser of a Product who is an end-user of such Product that the product is free from defects in material and workmanship, ordinary wear and tear excepted, for one year from date of purchase by such end-user. Further, Supplier warrants that the Product conforms to the NIOSH standards in effect on that date of sale. Supplier's sole obligation, and each end-users sole and exclusive remedy under this warranty, shall be limited to replacement of any Product covered by this warranty which proves in violation of such warranty during the period of effectiveness thereof. The warranty shall not apply to any Product which has been altered, misused, subjected to an accident, used contrary to Supplier's specifications and recommendations or used beyond its useful life or to any Product whose packaging has been altered or opened prior to purchase by the end-user. Supplier makes no representations or warranties concerning the application or effect of any international, federal, state or local laws, ordinances, rules or regulations which may now or hereafter pertain to the use or operation of the Product. EXCEPT AS EXPRESSLY PROVIDED IN THIS PARAGRAPH, SUPPLIER MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND, NATURE OR DESCRIPTION, EXPRESSED OR IMPLIED, INCLUDING, WITHOUT LIMITATION, ANY WARRANTY OF MERCHANTABILITY, FITNESS OF ANY OF THE PRODUCTS FOR ANY PARTICULAR PURPOSE, OR NON-INFRINGEMENT, AND HEREBY DISCLAIMS THE SAME.

In order to assert your rights under this warranty, please return the Product to your vendor together with a written description of the alleged defect and proof of purchase. Your distributor will forward the Product to Supplier for analysis of the validity of your claim.

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